

→ Compact and MR safe

MediGlasses allow participants with reduced visual acuity to see clearly within the scanner room and bore. These all-plastic glasses are suitable for use in MRI, fMRI, MEG, PET, CT and EEG environments. The slim, contour style of the adult-sized frame is particularly suitable for use with the latest multi-channel head coil designs, where there is often very limited space between the participant's face and the top of the coil.



→ Comfortable and easy to use

No previous optical knowledge is necessary for clear, quick, and safe vision correction. Lenses are easily exchanged with a unique thumb screw system. MediGlasses frames feature an adjustable head strap and face/nose padding for supreme comfort.



Also included in every MediGlasses set is a vision test chart, lens cleaning liquid and micro-fibre cleaning cloth, all packaged in a compact case to facilitate storage and handling.

→ Available in three versions

Each MediGlasses set offers a wide range of optical powers:

- **Full Set:** Twenty-six pairs of -6 to +6 dioptre lenses in 0.5 dioptre increments.
- **Standard Set:** Fifteen pairs of -5 to +3 dioptre lenses in 0.5 dioptre increments.
- **Mini Set:** Nine pairs of -5 to +3 dioptre lenses in 1.0 dioptre increments.



Cambridge Research Systems

Tel: +44 1634 720707

USA/Canada Toll Free: 1 866 846 2929

Email: enquiries@crsltd.com

www.crsltd.com

For more details:

www.crsltd.com/mediglasses



CAMBRIDGE RESEARCH SYSTEMS

M0111
M0105
M0106



→ Improve scan results

NoMoCo memory foam pillow support system products help patients complete brain imaging studies in research centres, hospitals or clinical settings with maximum comfort and with minimal movement. This results in the collection of high quality image output data that eliminates the need for repeat tests.

→ Designed for support

NoMoCo pillows are precision cut and hand made from the highest grade memory foam. Each twelve-piece set includes:

- **Head rest pads:** Placed in the cradle of the head coil for comfort.
- **Neck rolls:** Offer cervical support.
- **Wedge slide pads:** Offer lateral support and can also be used to hold headphones in place.
- **Square pads:** Placed around the top of the head to add further support, particularly for smaller subjects.
- **Head strap:** Restrains the forehead and provides sensory feedback regarding head placement and positioning.



→ Clinical and research coatings

NoMoCo pillows are offered with two types of coatings. The rugged *Healthcare* coating withstands repeated cleaning using hospital grade disinfectants.

The *Vinyl* coating is best suited for non-clinic research on healthy subjects. It's rated by ASTM standards to provide very good chemical resistance to various fluids, including Betadine, blood, acid, urea, alcohol and all-purpose cleaning products.



Cambridge Research Systems

Tel: +44 1634 720707

USA/Canada Toll Free: 1 866 846 2929

Email: enquiries@crsltd.com

www.crsltd.com

For more details:

www.crsltd.com/nomoco



CAMBRIDGE RESEARCH SYSTEMS

N1702
N1703